



# Upsee

Upsee can open up your child's world. It's a revelation for many kids with neuromuscular disorders. Some have been able to stand and hug their little brother or sister for the first time. Some have waved at neighbours on their first ever walk along their own street. And some crossed the playground hand-in-hand with their best friend for the very first time.



## THE MOBILITY DEVICE CHANGING PEOPLES LIVES WORLDWIDE

Upsee is an upright mobility device designed to allow children with motor impairment to stand and move with the help of an adult.

## A MOTHERS INVENTION

Debby's son, Rotem, was born with cerebral palsy. When he was just two-years-old, his therapists told Debby that her son didn't know what his legs were for. But neither Rotem or Debby were easily deterred. Debby began looking for ways to help Rotem stay active and use his arms and legs as much as possible. "I started finding ways of attaching him to me. I put straps around our legs. I tried different connections and materials until I got to the earliest version of the Upsee."

*They are just little steps – but they can mean the world to someone.*

## PARTICIPATION & INCLUSION

Upsee makes play and participation possible right now. It creates the possibility of social and emotional development by letting children interact with their peers at eye-level.

## COLOUR YOUR UPSEE

Your child has unique personality with unique abilities, and their Upsee Seat reflects that. Choose from two different colours in each size.

### X SMALL COLOURS



BLUE



PURPLE

### OTHER SIZES COLOURS



BLUE



PINK

## MOTIVATION & DEVELOPMENT

When your child is motivated to explore, move and weight-bear, it creates a real possibility for long-term development.

## OPPORTUNITIES FOR USE

Leading therapists endorse Upsee because it creates so many possibilities. It can be used for fun and play around the home or with your child's therapist during therapy sessions.

## WORK ON IMPORTANT SKILLS

Upsee can be used to work on a variety of mobility and ability skills, from weight-bearing and walking to improving head control.

## THE KIT CONTAINS

Upsee is made up of three components and the patented design uses durable, machine-washable and breathable fabrics so you and your child find it comfortable, reliable and easy to clean. These components include, the adult hip belt for the parent or carer to strap around their hips, the harness which holds the child securely in place, and the shared sandals.

## PRODUCT OVERVIEW

- Suitable for children aged 1-8 with a maximum user weight of 25kg/55lbs.
- The Upsee kit is made up of three parts: an adult hip-belt, a child harness and shared sandals.
- Facilitates gradual weight-bearing, standing and walking with the assistance of an adult.
- It is a revelation for many kids with neuromuscular disorders.



ADULT BELT



SANDALS



XS (BABY)  
HARNESS



CHILD  
HARNESS

Size	Shoulder to hip	Chest circumference	Lower hip circumference	Max user weight
XS	28-32cm	50-57cm	50-57cm	15kg
S	31-36cm	52-60cm	52-60cm	15kg
M	35-39cm	57-65cm	58-67cm	20kg
L	38-42cm	63-71cm	64-73cm	25kg